

# Molecular biology products for research and diagnostics



[bioron.de](http://bioron.de)

2009  
2010



## THERMOPHILIC DNA POLYMERASES

<b>Features of Thermophilic Polymerases</b>	<b>8</b>
<b>Selection table for Thermophilic Polymerases</b>	<b>9</b>
DFS-Taq DNA Polymerase	10
DFS-Blue Taq DNA Polymerase	11
Top-Taq DNA Polymerase	12
SuperHoTTaq	12
AntiTaq	13
DF-Pfu DNA Polymerase	14
DF-Psp DNA Polymerase	14
Klein-Taq DNA Polymerase	15
SNPase	16
Encyclo PCR Kit	17
Tersus PCR Kit	17
Tth DNA Polymerase	18

## MASTER MIXES

<b>Master Mixes</b>	<b>19</b>
<b>Features of Master Mixes</b>	<b>20</b>
<b>Selection table for Master Mixes</b>	<b>20</b>
DFS-Master Mix	21
<b>Hot Start PCR Master Mixes</b>	<b>21-23</b>
SuperHot Master Mix	21
HotRox Master Mix	22
dU-Hot Master Mix	23
<b>Uracil-DNA Glycosylase (UDG; related product)</b>	<b>23</b>
<b>Lyophilized Hot Start PCR Master Mixes</b>	<b>24-25</b>
Dry Master Mix blue	24
Dry Master Mix qPCR	25

## REVERSE TRANSCRIPTION

<b>Reverse</b>	<b>26</b>
<b>AMV reverse transcriptase</b>	<b>27</b>
<b>1 Tube RT-PCR Kit for qPCR</b>	<b>27</b>
<b>Plastic capillaries for Light Cycler</b>	<b>28</b>



<b>Diagnostic Kits for real-time PCR with CE-Conformity declaration</b>	<b>30-33</b>
Mycoplasma pneumoniae test system	30
MRSA screen - test system	31
Legionella test systems	32
Treponema pallidum (Tp) test system	33
<b>PCR Quantification Standards</b>	<b>34</b>
<b>S-Taq DNA Polymerase (qPCR) for test systems</b>	<b>34</b>
<b>GenePack DNA/RNA PCR tests systems (for R&amp;D only)</b>	<b>35-38</b>
GenePack DNA PCR tests	36
GenePack RNA PCR tests	38

#### DNA/RNA PURIFICATION

DNA/RNA purification with Ron's Line	40
Plasmid DNA purification Kits	41
Genomic DNA purification with Ron's Kits	42
Blood Genomic DNA Mini Kit	42
Tissue Genomic DNA Mini Kit	43
Fungal Genomic DNA Mini Kit	43
PCR-Pure DNA Kit	44
Gel Extraction Kit	44
Total RNA Kit	45
FastTri Reagent	45
OTP (one-tube-purification) Reagent	46
PCR-ready-Reagent (DNA Purification from blood samples)	46

#### PROTEIN PURIFICATION

Bio Ni-IDA Mini Kit	47
Bio Glutathione A for purification of GST-tagged Proteins	48
Total Protein Extraction Reagent for Mammalian Cells	49
Cytoplasmic and Nuclear Protein Extraction Kit	50



## 4 DNA AND DNA MARKERS

51-60

DNA Molecular Weight Markers Comparison Chart	52
DNA Ladder	53-55
50 bp DNA Ladder	53
100 bp DNA Ladder	53
100 bp extended DNA Ladder	54
1 kb DNA Ladder	54
Phage Lambda DNA/BstE II Marker	55
Phage Lambda DNA/Hind III Marker	56
Phage Lambda DNA/Sty I Marker	56
pUC 19 plasmid	57
pUC 19/Msp I DNA Marker	57
Loading Buffer II (6x)	58
Loading Buffer IV (10x)	58
Phage Lambda DNA	59
Phage T7 DNA	59
Human DNA	60
Cot I Human DNA	60

## 5 DEOXYNUCLEOSIDE 5'-TRIPHOSPHATES

61-72

### STANDARD DNTPS

2'-Deoxynucleoside 5'-triphosphates	62
dNTP Mix	63
Set of 4 dNTPs (dATP, dCTP, dGTP, dTTP)	63
dATP	64
dCTP	64
dGTP	65
dTTP	65
dUTP	66

### MODIFIED DNTPS

Bio-11-dUTP	66
Deoxynucleoside triphosphates with aminolinkers	67
Aminohexyl-dCTP	67
Aminobutyl-dATP	68
Aminoallyl-dUTP	68
Flu-12-dUTP	69
TAMRA-12-dUTP	69

### CHAIN TERMINATORS

ddATP	70
ddCTP	70
ddGTP	71
ddUTP	71
Cy5-11-ddUTP	72



Restriction endonucleases _____	74
Alphabetic list _____	75
List of Prototypes _____	76

List of Bioron's endonucleases _____	77-88
--------------------------------------	-------

Alu I	AG↓CT
ApaI	G↓TGCAC
Asu II	TT↓CGAA
BamH I	G↓GATCC
Bcl I	T↓GATCA
Bgl I	GCCNNNN↓NGGC
Bgl II	A↓GATCT
BseA I	T↓CCGGA
BseB I	CC↓(AT)GG
BseC I	TA↓CGAT
BshF I	GG↓CC
BsiS I	C↓CGG
BssA I	R↓CCGGY
BstE II	G↓GTNACC
CspA I	A↓CCGGTT
EcoR I	G↓AATTC
EcoR V	GAT↓ATC
Hind III	A↓AGCTT
Hinf I	G↓ANTC
Hpa I	GTT↓AAC
Kpn I	GGTAC↓C
Mbo I	↓GATC
MspC I	C↓TTAAG
Nae I	GCC↓GGC
Nco I	C↓ATGG
Nhe I	G↓CTAGC
Not I	GC↓GGCCGC
Nru I	TCG↓CGA
Pst I	CTGCA↓G
Pvu II	CAG↓CTG
Rsa I	GT↓AC
Sal I	G↓TCGAC
Sau3A I	↓GATC
Sca I	AGT↓ACT
Sfi I	GGCCNNNN↓NGGCC
SgrB I	CCGC↓GG
Sla I	C↓TCGAG
Sma I	CCC↓GGG
SnaB I	TAC↓GTA
Sph I	GCATG↓C
SseB I	AGG↓CCT
Ssp I	AAT↓ATT
SstI	GAGCT↓C
Sty I	C↓CWWGG
Taq I	T↓CGA
Xba I	T↓CTAGA

Restriction endonucleases activity chart _____	89
Buffers _____	90



**7 ENZYMES AND FINE CHEMICALS**

91-98

**ENZYMES**

Proteinase K	92
T4 DNA Ligase	92
T4 RNA Ligase	93
T4 Polynucleotide Kinase	93
Exonuclease III	94
Tli Inorganic Pyrophosphatase	94
DNA Polymerase I, Large (Klenow) Fragment	95
φ 29 DNA Polymerase	95

**FINE CHEMICALS**

Exo Resistant Random primer	96
Oligo(dT) <sub>15</sub> Primer	96
Random Primers	96
Linear Polyacrylamide Solution	97
Agarose Molecular Biology Grade	97
X-Gal	98
IPTG	98

**8 DEVICES**

99-104

Introduction	100
Gene Machine	101
DNA-diagnostics with Gene Machine	101
SpeedScan	102
Spectrophotometer A-007	104

**9 CONSUMABLES**

105-114

**LABORATORY PLASTIC**

Polycarbonate Capillaries for LightCycler	106
Universal tips	107
Filtertips	108
Test tubes	109
Thin-walled PCR-tubes	109

**DIALYSIS**

High Quality Regenerated Cellulose Dialysis Membranes	110-114
Cellu•Sep T-Series Dialysis Membranes	110
Cellu•Sep H1-Series Dialysis Membranes	112
QuixSep	114



## BIORON Germany

<b>Tel.</b>	+49 621 5720 915
<b>FAX</b>	+49 621 5720 916
<b>E-mail</b>	info@bioron.net
<b>Post Bioron GmbH</b>	Rheingönheimer Str. 36 D-67065 Ludwigshafen Germany

## BIORON Singapore

<b>Tel.</b>	+65 6896 8063
<b>FAX</b>	+65 6896 8092
<b>E-mail</b>	sales@bioron.net
<b>Post Bioron PVT. LTD.</b>	Room 9 No. 51 Bukit Batok Crescent Unity Centre #07-05 Singapore 658077

## BIORON Poland

<b>Tel.</b>	+48 426 7704 57
<b>FAX</b>	+48 426 7704 58
<b>E-mail</b>	Grzegorzczuk@bioron.net
<b>Post Bioron Sp. z.o.o.</b>	Puszkina str. 80 92-516 Lodz Poland

**To facilitate your order, please provide the following information**

- Institution or company customer number
- Name, telephone number
- Billing address
- Shipping address
- Valid Purchase Order number
- Catalogue number, description and quantity of the products

**Place your order online at:  
24 hours-ONLINE-WEB-SHOP**

[www.bioron.net/shop](http://www.bioron.net/shop)

BIORON GmbH  
Tel.: +49 621 5720915  
Fax: +49 621 5720916  
E-mail: [info@bioron.net](mailto:info@bioron.net)

BIORON PTE LTD (Singapore)  
Tel.: +65 6896 8063  
Fax: +65 6896 8092  
E-mail: [sales@bioron.net](mailto:sales@bioron.net)

BIORON Sp. z.o.o. (Poland)  
Tel.: +48 426 7704 57  
Fax.: +48 426 7704 58  
E-mail: [grzegorzcyk@bioron.net](mailto:grzegorzcyk@bioron.net)

[www.bioron.de](http://www.bioron.de)