

# HMC CYCLONE

## Magnetic Stirrers

- Ⓢ Electronic magnetic stirrer without moving parts
- Ⓢ Remote controllable single-point stirrer in aluminium housing
- Ⓢ Hermetically sealed
- Ⓢ For using in devices or in water baths
- Ⓢ Rotating speed 100-1000 rpm
- Ⓢ Integrated M4 threaded bore holes for attaching stirrer in analysis devices tripods and so on
- Ⓢ Indestructible



HMC EUROPE GmbH  
Labor- und Sterilisationstechnik  
Hafing 21  
D - 84549 Engelsberg - Germany

Tel. +49 (8634) 625 994 /-5  
Fax +49 (8634) 625 996  
info@hmc-europe.com

CYCLONE 1-60 - footprint 60 mm Order No. 90003  
CYCLONE 1-100 - footprint 100 mm Order No. 90004  
CYCLONE 1-150 - footprint 150 mm Order No. 90005  
incl. CYCLOPOWER control unit

[www.hmc-europe.com](http://www.hmc-europe.com)