

Labotron

The handy bench-scale orbital shaker



- ▶ Flexible
- ▶ Reliable
- ▶ Robust
- ▶ Corrosion resistant



Reliable shaking of small loads

► Flexible

Labotron is at home on the laboratory bench or in an incubator.

► Reliable

The sealed-for-life bearing lubricated with a high-temperature grease ensures a long and productive working life whatever the environment.

► Robust

The unit has a strong housing made from composite material. The shaking table is made from a solid plastic. This makes the whole easy to clean and resistant to abrasion and a wide range of chemicals.

► Corrosion resistant

The composite housing is corrosion resistant and ensures easy cleaning.



Loading Capacity

| | Size of flask | Quantity |
|--|---------------|----------|
| For operation using accessories without the need for the special screws, a table can be fitted. This is of type ATW 300 X 370 mm and can be universal, fixed or with adhesive matting. | 50 | 35 |
| | 100 | 23 |
| | 250 | 12 |
| | 500 | 8 |

Technical Data

| | |
|----------------------|-------------------------|
| Shaker Speed | 40–400 RPM |
| Throw | ø 12,5 mm |
| Permissible | |
| Ambient T° | 4°C up to 60°C |
| Ambient humidity | up to max. of 80% rF |
| Weight without table | Approx. 9 kg |
| Dimensions W x D x H | 39 x 32 x 9 cm |
| Power supply | 115/230V ± 10%, 50/60Hz |

Infors AG
Headoffice, Switzerland

Rittergasse 27
CH-4103 Bottmingen
T +41 (0)61 425 77 00
F +41 (0)61 425 77 01
headoffice@infors-ht.com

Infors UK Ltd

The Courtyard Business
Centre, Dovers Farm
Lonesome Lane, Reigate
Surrey, UK-RH2 7QT
T +44 (0)1737 22 31 00
F +44 (0)1737 24 72 13
infors.uk@infors-ht.com

Infors Sarl

6, rue Marcel Paul
FR-91742 Massy Cedex
T +33 (0)169 30 95 04
F +33 (0)169 30 95 05
infors.fr@infors-ht.com

Infors Italia S.r.l.

Via Barbara Melzi, 13
I-20025 Legnano (MI)
T +39 0331 59 47 67
F +39 0331 54 40 99
infors.it@infors-ht.com

Infors GmbH

Tannenweg 2
D-85254 Einsbach
T +49 (0)8135/83 33
F +49 (0)8135/83 20
infors.de@infors-ht.com

Infors Benelux BV

Logistiekweg 20-02
7007 CJ Doetinchem
P.O. Box 613
NL-7000 AP Doetinchem
T + 31(0)314-36 44 50
F + 31(0)314-37 80 76
infors.bnl@infors-ht.com

For more information please visit:
www.infors-ht.com