

# Silver Staining Solution


A rapid and highly sensitive silver stain system for protein detection

## ■ Features

- A rapid and easy method for staining proteins in polyacrylamide gels
- 30 times more sensitive than Coomassie staining method
- No background and free from dark spots
- Complete destaining of gels within 15 minutes
- Unique buffers and protocol for preparation of protein sample for MALDI-TOF-MS analysis (no glutaraldehyde included)
- Sensitivity: Detects 0.1 ng of BSA
- No. of Applications: 20 mini-gels
- Storage : 6 months at room temperature

## ■ Kit components

- Sensitizer ..... 50 ml
- Stainer A ..... 5 ml
- Stainer B ..... 5 ml
- Developer A ..... 50 ml
- Developer B ..... 5 ml



Cat.No.	Description	Size
K14040D	EzWay™ Protein-Silver Staining Kit	20 mini-gels