

DNA/RNA Purification Kit


Fast and efficient isolation & extraction of DNA or RNA

■ Features

- Easy-to-Use protocol for isolation of DNA or RNA from various samples
- Silica-based membrane based technology in the spin column
- High yield and purity of isolated DNA or RNA
- Time saving procedure



■ Principle



DNA or RNA released from sample treated by lysis buffer is bound to a silica gel membrane in the presence of high concentration of chaotropic salt. The purified genomic DNA on the membrane is eluted in the elution buffer.

Cat.No.	Description	Size
K11738	EzWay™ Genomic DNA Kit, Virus	50T
K11739	EzWay™ Genomic DNA Kit, Virus	200T
K11740	EzWay™ Genomic DNA Kit, Blood/Cultured Cells	50T
K11741	EzWay™ Genomic DNA Kit, Blood/Cultured Cells	200T
K11742	EzWay™ Genomic DNA Kit, Plant	50T
K11743	EzWay™ Genomic DNA Kit, Plant	200T
K11730	EzWay™ Genomic DNA Kit, Multi (Animal/Bacteria/Yeast)	50T
K11731	EzWay™ Genomic DNA Kit, Multi (Animal/Bacteria/Yeast)	200T
K33305	EzWay™ Buccal Swab DNA Isolation Kit	50T
K33310	EzWay™ Buccal Swab DNA Isolation Kit	100T
K11701	EzWay™ Plasmid DNA Kit	50T
K11702	EzWay™ Plasmid DNA Kit	200T
K11760	EzWay™ PCR Clean-Up Kit	50T
K11761	EzWay™ PCR Clean-Up Kit	200T
K11770	EzWay™ Gel Extraction Kit	50T
K11771	EzWay™ Gel Extraction Kit	200T
K11755	EzWay™ miRNA Isolation Kit	100T
K11750	EzWay™ Total RNA Isolation Kit	50T
K11751	EzWay™ Total RNA Isolation Kit	200T