

[Request Quote](#) | [Specifications](#) | [Order Information](#)

Centrifuges :: Minicentrifuges

Spectrafuge Mini Centrifuge

Join the mini revolution



- Supplied with both standard microtube and strip tube rotors
- Ideal for quick spin downs and microfiltration
- Starts and stops in seconds
- Compact design



Product Information

This personal sized centrifuge is ideal for microfiltration and quick spin downs from the walls and caps of microcentrifuge tubes



Requiring less than six inches of bench space, the Spectrafuge Mini is supplied with rotors and adapters to accommodate 1.5 ml, 0.5 ml and 0.4 ml tubes as well as 0.2 ml strips and individual tubes. The standard and strip tube rotors are easily interchanged.

The easy-open, translucent lid utilizes a stainless steel hinge pin for durability. Five different lid colors are available including gray, blue, teal, purple and red.

An on/off switch is located on the side of the Spectrafuge Mini for starting and stopping operation. Alternatively, with the switch in the on position, the centrifuge can be started and stopped by closing and opening the lid.



Back to Top

Product Specifications

Maximum speed/RCF	6000rpm/2000 x g
Maximum capacity	
Strip tube rotor	2 x 0.2 ml strips
Standard tube rotor	6 x 1.5/2.0 ml
Dimensions (W x D x H)	5 x 5 x 4.4 inches/12.7 x 12.7 x 11.2 cm
Weight	1 lb/0.45 kg
Electrical*	120V~, 60 Hz
Warranty	12 Months
1 Average max rcf in the strip tube rotor is 1177 x g *Other voltages available upon request	

Back to Top

Ordering Information

Catalog Number	Item Description	Price
C1301	Spectrafuge Mini complete with 1.5/2.0 ml rotor, strip tube rotor, 0.5 ml and 0.4 ml adapters, smoke gray lid	Request Quote
C1301-B	Spectrafuge Mini, as above, blue lid	Request Quote

C1301-T	Spectrafuge Mini, as above, teal lid	Request Quote
C1301-P	Spectrafuge Mini, as above, purple lid	Request Quote
C1301-R	Spectrafuge Mini, as above, red lid	Request Quote
C1300-RT	Replacement rotor for 6 x 1.5/2.0 ml tubes	Request Quote
C1300-RTS	Replacement rotor for 0.2 ml strip tubes or individual 0.2 ml tubes	Request Quote
C1205	Individual adapters for 0.5/0.6 ml tubes, pk of 6	Request Quote
C1206	Individual adapters for 0.4/0.25 ml tubes, pk of 6	Request Quote
C1222	Individual adapters for 0.2 ml thermal cycling tubes, pk of 6	Request Quote

*When ordering 220V units, please add -220V to the catalog number. 220V units supplied with European plug ends.

[Back to Top](#)

Download related documents

[Spectrafuge Mini Centrifuge Instruction Manual](#)

Labnet International - International Products | Ph: +1-732-417-0700 | Fax: +1-732-417-1750 | international@labnetlink.com

[Products](#)

[What's New](#)

[Contact Us](#)

[Register
Product](#)

[Ordering](#)

[Search](#)

[OEM
Offerings](#)

[Service](#)

[Labnet Main](#)