



## **Smart / Pink Minis 300V Power Supply**

# Smart / Pink Minis 300V Power Supply

The Minis 300V power supplies are microprocessor-controlled power supplies with universal rated voltage with up to 300 voltage output. The ultracompact sized unit provides power needed to run horizontal agarose or vertical polyacrylamide electrophoresis. The Minis 300V series provide 6 different constant voltage selections in timed or continuous operation. Easy to operate, dual output terminators, 6 constant voltages, and compact size are the features. The Minis 300V series is also designed and constructed to the most rigorous safety standards. The packages with 10 x 10 cm vertical electrophoresis or mild sized horizontal electrophoresis systems are valuable choices for educational market and personal use.

## Features

- >50V or 100V or 150V or 200V or 250V or 300V constant voltage output
- >400mA maximum current
- >60W maximum power
- >Two pairs of outlet terminator
- >Timer with alarm function
- >Safety device
- >Compact size and lightweight
- >Competitive pricing choice
- >Packages with mild sized horizontal electrophoresis systems and vertical electrophoresis systems



CE

Cat. No. :Minis 300VS



CE

Cat. No. :Minis 300VP

Specification	Minis 300VS	Minis 300VP
Appearance	Blue color faceplate	Pink color faceplate
Max. Voltage	300V	
Max. Current	400mA	
Max. Watt	60W	
Output Type	Constant voltage	
Control	Microprocessor controller	
Voltage Selection	50V, 100V, 150V, 200V, 250V, 300V	
Terminator Pairs	2 Pairs	
Timer	1~999 mins with alarm, continuous	
Safety Device	No load detection; shrouded plugs and sockets	
Operation Temperature	Ambient to 40°C	
Unit Dimension	115 x 150 x 80mm (W x L x H)	
Construction Material	Frame retard ABS and aluminum rare and bottom plates	
Weight	approx. 0.6kg	
Rated Voltage	100~240V	
Application	Sufficient for mild sized horizontal electrophoresis systems and vertical electrophoresis systems	

### packages available



ME15-15



MV-10DSYS

ME15-300VS



ME15-300VP




MV10-300VS




MV10-300VP





Two pairs of outlet terminator



Compact size and lightweight

Universal rated voltage

100~240V

**Ordering information**

Minis 300VS Smart Minis 300V Power Supply

Minis 300VP Pink Minis 300V Power Supply

ME15-300VS Package of ME15-15 and Smart Minis 300V

MV10-300VS Package of MV-10DSYS and Smart Minis 300V

ME15-300VP Package of ME15-15 and Pink Minis 300V

MV10-300VP Package of MV-10DSYS and Pink Minis 300V



## Contact Information

[www.majorsci.com](http://www.majorsci.com)

[info@majorsci.com](mailto:info@majorsci.com)

19959 Sea Gull Way  
Saratoga, CA 95070  
U.S.A

No.7, Hwai-Der Street,  
Pan-Chiao City, Taipei Hsien 22045,  
Taiwan

T/1-408-366-9866  
F/1-408-446-1107

T/886-2-2256-6980  
F/886-2-2256-6942

*Innovative Life Sciences Tools*