



PERLONG GROUP
Perlong Medical Equipment
Co., Ltd.

Perlong Group-Cherish your life, Cherish your health!

Clinical Lab Devices

PUS-2018G Biochemistry Analyzer



★Specifications★

- ◇ Wavelength Range: 340nm~750nm
- ◇ Measurement Range: 0.000~3.000Abs
- ◇ Repeatability: ≤ 0.005
- ◇ Linearity: $\leq 1.0\%$
- ◇ Cross-contamination: $\leq 1.0\%$
- ◇ Temperature Accuracy: $\pm 0.2^{\circ}\text{C}$
- ◇ Colorimetric Pool: Titanium alloy
- ◇ Test Method: Absorbency, End-poit, Kinetic, Two-point, Standard Curve, Control Tube, Facto

★Features★

- ◇ 7.5" colorful touch screen, can be connected to a mouse, Windows operating interface
- ◇ Can store 200~500 test results
- ◇ Can show real-time reaction curve, and print calibration curve, kinetic reaction curve, quality-control curve and general report
- ◇ While not working, the instrument will automatically enter half-voltage sleep status to save energy

PERLONG GROUP

Add: No.28-1 Dajiaochang Road, Nanjing, 210022 China

Tel: 86-25-52635350 52635360 Fax: 86-25-52635356

<http://www.perlong.com> E-mail: nanjingpt@perlong.com.cn