

INCUDRIVE R

Roller system

- With drive unit and roller frame, without roller drums
- To slowly rotate drums with test-tubes and bottles
- Can be combined with INCUDRIVE H

Rotational speed of drums 0.1 to 4.8 rpm as well as inclination of frame are continuously adjustable.

For use in an incubating room or any other temperature-controlled environment.

Drums for test-tubes with 100 and 160 mm length available. Custom-made roller drums for other vessels can be manufactured on request.



INCUDRIVE R
and drums (optional)

Technical Data

Dimensions outer (wxhxd):	650 x 260 x 530 mm (without drum)
Power:	230 V, 50/60 Hz, 80 W
Weight:	approx. 15 kg

Ordering Information

Type	Cat.-No.
INCUDRIVE R Roller-System (Drive unit with roller frame, without drum)	3.064 001
Accessories	
Drum R for 7 roller bottles Ø 120 mm (max. 1 each per roller frame)	3.056 071
drum K10 for 210 test-tubes 100x16 mm (max. 3 each per roller frame)	3.056 101
Drum K16 for 210 test-tubes 160x16 mm (max. 2 each per roller frame)	3.056 161
drum M for other vessels (custom-made) (max. 2 each per roller frame)	3.056 311