

# INCUDRIVE H

## Incubating hood

- with air-circulation system
- bottom open
- can be combined with INCUDRIVE R

For use with various kinds of experiments or apparatus on even surfaces.

Temperature range from +5 °C above room temperature to +50 °C, temperature continuously adjustable and stabilized, forced air-circulation.

Transparent walls permit permanent monitoring of experiment.

Soft rubber lining at the bottom for passage of power cords etc.



INCUDRIVE H

### Technical Data

Dimensions inner (wxhxd):	600 x 585 x 480 mm
Dimensions outer (wxhxd):	650 x 750 x 570 mm
Power:	230 V, 50/60 Hz, 280 W
Weight:	approx. 10,5 kg

### Ordering Information

Type	Cat.-No.
INCUDRIVE H Incubating Hood	3.065 001