

# ALU-STECK-caps

**Purealuminium.** Matt finished. Sterilizable

For placing loosely on top of glass vessels.  
Sterile but non-hermetic closing of vessels.  
Without spring.

Can also be used as open sample containers.



ALU-STECK-caps

## Ordering Information

Inner-Ø mm	Height mm	Qty.	Cat.-No.
11	25	100	<b>3.620 113</b>
12	25	100	<b>3.620 213</b>
13	25	100	<b>3.620 313</b>
14	25	100	<b>3.620 413</b>
15	25	100	<b>3.620 513</b>
16	30	100	<b>3.620 613</b>
17	30	100	<b>3.620 713</b>
18	30	100	<b>3.620 813</b>
20	30	100	<b>3.620 913</b>
22	30	100	<b>3.621 113</b>
28	30	100	<b>3.621 213</b>
32	30	100	<b>3.621 313</b>
39	30	10	<b>3.621 413</b>
44	40	10	<b>3.621 513</b>
48	40	10	<b>3.621 613</b>
55	40	10	<b>3.621 713</b>
58	40	10	<b>3.621 813</b>