

CAP-O-TEST-caps

Aluminium. Anodized in color. Anodic oxidised. Sterilizable.

Sterile but non-hermetic closing of vessels.

Springs made of chromium-nickel'-steel, which clamp against the inner wall of the vessel.

Special colors available on request.



CAP-O-TEST-caps

Ordering Information

Vessel with \varnothing mm	Qty.	Cat.-No. silver	Cat.-No. blue	Cat.-No. red	Cat.-No. yellow	Cat.-No. green	Cat.-No. brown	Cat.-No. black	Cat.-No. violet
12/13	100	3.626 113	3.626 123	3.626 133	-	-	-	-	-
15/16	100	3.626 313	3.626 323	3.626 333	3.626 343	3.626 353	3.626 363	3.626 373	3.626 383
16/17	100	3.626 413	3.626 423	3.626 433	3.626 443	3.626 453	3.626 463	3.626 473	3.626 483
17/18	100	3.626 513	3.626 523	3.626 533	-	-	-	-	-
19/20	100	3.626 613	3.626 623	3.626 633	-	-	-	-	-
24/26	50	3.626 813	3.626 823	3.626 833	-	-	-	-	-