

LABOCAP-caps

with/without handle

Aluminium. Anodized in color. Anodic oxidised. Sterilizable.

For sterile but non-hermetic closing of vessels. Springs made of chromium-nickel steel, which clamp against the outer wall of the vessel.

Other sizes and special colors available on request.



LABOCAP-caps with handle

Ordering Information

(with handle)

| Vessel with \varnothing mm | Qty. | Cat.-No. silver | Cat.-No. blue | Cat.-No. red |
|------------------------------|------|-----------------|---------------|--------------|
| 12/13 | 100 | 3.623 113 | 3.623 123 | 3.623 133 |
| 14/15 | 100 | 3.623 213 | - | - |
| 15/16 | 100 | 3.623 313 | 3.623 323 | 3.623 333 |
| 17/18 | 100 | 3.623 413 | 3.623 423 | 3.623 433 |
| 19/20 | 100 | 3.623 513 | 3.623 523 | 3.623 533 |



LABOCAP-caps without handle

(without handle)

| Vessel with \varnothing mm | Qty. | Cat.-No. silver | Cat.-No. blue | Cat.-No. red | Cat.-No. yellow | Cat.-No. green | Cat.-No. black |
|------------------------------|------|-----------------|---------------|--------------|-----------------|----------------|----------------|
| 09/10 | 100 | 3.624 113 | - | - | - | - | - |
| 12/13 | 100 | 3.624 213 | 3.624 223 | 3.624 233 | 3.624 243 | 3.624 253 | - |
| 14/15 | 100 | 3.624 313 | 3.624 323 | 3.624 333 | - | - | - |
| 15/16 | 100 | 3.624 413 | 3.624 423 | 3.624 433 | 3.624 443 | 3.624 453 | 3.624 473 |
| 17/18 | 100 | 3.624 513 | 3.624 523 | 3.624 533 | 3.624 543 | 3.624 553 | - |
| 19/20 | 100 | 3.624 613 | 3.624 623 | 3.624 633 | - | - | - |
| 21/23 | 100 | 3.624 713 | - | - | - | - | - |
| 24/26 | 100 | 3.624 813 | - | - | - | - | - |
| 28/30 | 100 | 3.624 913 | - | - | - | - | - |