

ERLENMEYER-CAP-caps

Aluminium. Anodized in color. Anodic oxidised. Sterilizable.

Sterile but non-hermetic closing of vessels.

Springs made of chromium-nickel-steel, which clamp against the outer wall of the vessel.

Special caps for culture flasks and culture media bottles with straight neck of \varnothing 38 mm.

Other sizes are available on request.



ERLENMEYER-CAP-caps

Ordering Information

Vessels with \varnothing mm	Qty.	Cat.-No. <u>silver</u>	Cat.-No. <u>blue</u>	Cat.-No. <u>red</u>
37/39	10	3.627 113	3.627 123	3.627 133