

Culture tubes

with thread and screw-cap

Soda lime-glass. Sterilizable.

For cultivation and for transporting, dispatch and long term storage.

Complete with aluminium screw-cap (autoclavable/sterilizable)
or plastic screw-cap (steam-autoclavable with 121 °C/sterilizable)
and rubber gasket.

Size 5 ml with frosted labeling space.



culture tubes with screw cap

Ordering Information

Size mm	Volume ml	Qty.	Cat.-No.
with aluminium screw-cap			
35x12	approx. 3	100	3.562 123
50x14 (with frosted labeling space)	approx. 5	100	3.562 143
with plastic screw-cap			
50x14 (with frosted labeling space)	approx. 5	100	3.563 143
45x27	approx. 15	100	3.563 423
50x27	approx. 20	100	3.563 523
75x27	approx. 30	100	3.563 623