

Magnetic stirrer ROTAMIX S-10 with stirring plate  $\varnothing$  135 mm is designed for stirring liquids in laboratories. ROTAMIX S-10 features a quiet motor drive with analog electronics. Speed is fully adjustable. Speed remains constant, even under varying loads. Stirring plate is made from stainless steel.

---



### MAIN FEATURES

Ideal for stirring applications  
Max. capacity 10 liters  
Stainless steel plate  
Convenient, small footprint design  
Analog speed

### TECHNICAL SPECIFICATIONS

Model	Rotamix S-10
Code No.	691.000
Voltage	230V/50Hz or 115V/60Hz
Power	40 VA
Speed range	200 – 1200 rpm
Stirring plate dimension	$\varnothing$ 135 mm
Max. capacity	10 liters
Ambient operating range	4° to 65°C
Dimensions (WxDxH)	170 x 265 x 135 mm
Net weight	4,2 kg

*We reserve the right to alter specification details etc. without prior notice or liability!*

---

### **TEHTNICA ŽELEZNIKI d.o.o.**

Podjetje precizne mehanike in elektronike, Na Plavžu 79, 4228 Železniki, Slovenia  
Tel.: +386 4 5117 500, Fax.: +386 4 5117 501, Internet: [www.tehtnica.si](http://www.tehtnica.si), E-mail: [info@tehtnica.si](mailto:info@tehtnica.si)

---