



Electric Furnace Technology

Tube Furnaces VTF choice of 1 - 2 or 4 tube models from 1200°/ 1600°C



Model	Max Temp	Max Cont. Working Temp.	Work Tube Bore (mm)	Heated Length (mm)	Extended Dimensions (mm)			Number of Tubes	Power (kW)	T/CPL Type
					Height	Width	Depth			
VTF1	1200°C	1150°C	25mm	150	560	400	340	1	1.5	K
VTF2	1200°C	1150°C	25mm	150	560	400	340	2	1.5	K
VTF3	1200°C	1150°C	25mm	150	560	400	340	4	1.5	K
VTF4	1400°C	1400°C	25mm	150	710	470	410	1	2.25	R
VTF5	1400°C	1400°C	25mm	150	710	470	470	2	2.25	R
VTF6	1400°C	1400°C	25mm	150	780	470	470	4	2.25	R
VTF7	1600°C	1600°C	25mm	150	550	520	460	1	3.5	R
VTF8	1600°C	1600°C	25mm	150	770	520	520	2	3.5	R
VTF9	1600°C	1600°C	25mm	150	830	520	520	4	3.5	R

* Please note larger bores available on request

Vecstar Ltd

Units 11 & 12 Foxwood Road Dunston Trading Estate Chesterfield S41 9RF
 Tel: 01246 260094 Fax: 01246 450213 Email: enquiries@vecstar.co.uk